

















BIOPLATE

CRANIAL FIXATION

1.5 MM SYSTEM




1.5MM PLATES

	Product Name	Description	Quantity	Code
	2-hole Straight Plate, Small	11 mm	1	81-3018
			5	81-2000
	2-hole Straight Plate, Large	14 mm	1	81-3131
			5	81-3135
	2-hole Straight Compression Plate	14 mm	1	81-3144
	4-hole Straight Plate	19 mm	1	81-2001
	4-Hole Boxcar Plate	10 x 10 mm	1	81-2007
	4-Hole Boxcar Plate, Large	14 x 14 mm	1	81-3136
	6-Hole Boxcar Plate	8 x 17 mm	1	81-2008
	8-Hole Boxcar Plate	10 x 24 mm	1	81-2009
	Y Plate	18 mm	1	81-2002
	X Plate	17 mm	1	81-2003
	T Plate	15 mm	1	81-2004
	Flexible Burr Hole Cover, Small	12 mm	1	81-2053
	Flexible Burr Hole Cover, Large	18 mm	1	81-2054
	Access Burr Hole Cover	18 mm	1	81-3180

Contact Bioplate to learn more about our complex and customizable plates not shown in this catalogue. 888-BIOPLATE


1.5MM SCREWS



Size	Quantity	Code
3 mm long	1	81-2065
	5	81-2071
4 mm long	1	81-2066
	5	81-2072
5 mm long	1	81-2067
	5	81-2073


1.5MM INSTRUMENTATION AND ACCESSORIES

SCREWDRIVER BLADE



Description	Quantity	Code
Screwdriver Blade	1	81-3101

SCREWDRIVER HANDLE



Description	Quantity	Code
Screwdriver Handle	1	81-3111

FORCEPS



Description	Quantity	Code
Forceps - Angled	1	81-2024

CONTAINER



Description	Quantity	Code
Plate & Screw Container	1	81-4205







1.5MM COMPACT PLUS STARTER KIT (81-3340)



Bioplate’s 1.5mm Compact Plus Starter Kit is a combination of the most frequently used plates, screws and instruments. It’s designed to be an introduction to our cranial fixation system. Delivered in our Compact Plus container, Bioplate provides a unique and flexible approach to rigid fixation. The Compact Plus container can be combined with other Bioplate containers for a stackable, compact footprint.

1.5MM STERILE KITS



Quantity Per Kit		Description	Code
3		2-hole Straight Plate, Large	81-3573 (BOX OF 10 KITS)
6		4 mm long Self-Drilling Screw	
Quantity Per Kit		Description	Code
1		X Plate	81-3578 (BOX OF 5 KITS)
6		4 mm long Self-Drilling Screw	
Quantity Per Kit		Description	Code
1		Flexible Burr Hole Cover, Large	81-3580 (BOX OF 5 KITS)
6		4 mm long Self-Drilling Screw	

Contact Bioplate to learn more about the contents of our **COMPACT PLUS KIT** and our customizable **STERILE KITS** not shown in this catalogue. 888-BIOPLATE

ORDER INFORMATION

888-BIOPATE

info@bioplate.com

PLATES - Ø1.5 mm System

81-3018	Straight Plate - 1.5 mm - 2 hole - 11 mm	Qty 1
81-2000	Straight Plate - 1.5 mm - 2 hole - 11 mm	Qty 5
81-3131	Straight Plate - 1.5 mm - 2 hole - 14 mm	Qty 1
81-3135	Straight Plate - 1.5 mm - 2 hole - 14 mm	Qty 5
81-3144	Straight Compression Plate - 1.5 mm - 2 hole - 14 mm	Qty 1
81-2001	Straight Plate - 1.5 mm - 4 hole - 19 mm	Qty 1
81-2007	Boxcar Plate - 1.5 mm - 4 hole - 10 x 10 mm	Qty 1
81-3136	Boxcar Plate - 1.5 mm - 4 hole - 14 x 14 mm	Qty 1
81-2008	Boxcar Plate - 1.5 mm - 6 hole - 8 x 17 mm	Qty 1
81-2009	Boxcar Plate - 1.5 mm - 8 hole - 10 x 24 mm	Qty 1
81-2002	Y Plate - 1.5 mm - 5 hole - 18 mm long	Qty 1
81-2003	X Plate - 1.5 mm - 6 hole - 17 mm long	Qty 1
81-2004	T Plate - 1.5 mm - 5 hole - 15 mm long	Qty 1
81-3180	Access Burr Hole Cover - 1.5 mm - 5 hole - 18 mm	Qty 1
81-2053	Flexible Contoured Burr Hole Cover - 1.5 mm - 12 mm	Qty 1
81-2054	Flexible Contoured Burr Hole Cover - 1.5 mm - 18 mm	Qty 1

STERILE KIT - Ø1.5 mm System

81-3573	Straight Plate - 1.5 mm - 2 hole - 14 mm 1.5 mm diam. Self-Drilling Screw - 4 mm long	(BOX OF 10 KITS)
81-3578	X Plate - 1.5 mm - 6 hole - 17 mm long 1.5 mm diam. Self-Drilling Screw - 4 mm long	(BOX OF 5 KITS)
81-3580	Flexible Contoured Burr Hole Cover - 1.5 mm - 18 mm 1.5 mm diam. Self-Drilling Screw - 4 mm long	(BOX OF 5 KITS)

SCREWS - Ø1.5 mm System

81-2065	1.5 mm diam. Self-Drilling Screw - 3 mm long	Qty 1
81-2066	1.5 mm diam. Self-Drilling Screw - 4 mm long	Qty 1
81-2067	1.5 mm diam. Self-Drilling Screw - 5 mm long	Qty 1
81-2071	1.5 mm diam. Self-Drilling Screw - 3 mm long	Qty 5
81-2072	1.5 mm diam. Self-Drilling Screw - 4 mm long	Qty 5
81-2073	1.5 mm diam. Self-Drilling Screw - 5 mm long	Qty 5

ACCESSORIES - Ø1.5 mm System

81-3101	Screwdriver Blade - 1.5 mm diam.
81-3111	Screwdriver - Handle only
81-2024	Forceps - Angled
81-4205	1.5 mm Plate & Screw Container

