



BIOPLATE

CRANIAL FIXATION

1.9 MM SYSTEM



1.9MM PLATES


	Product Name	Description	Qty	Code
	2-hole Straight Plate, Small	12 mm	1	81-3048
			5	81-3027
	2-hole Straight Plate, Large	15 mm	1	81-3020
			5	81-2037
	2-hole Straight Compression Plate	15 mm	1	81-3148
	4-hole Straight Plate	20 mm	1	81-3028
	4-Hole Boxcar Plate	11 x 11 mm	1	81-3034
	4-Hole Boxcar Plate, Large	15 x 15 mm	1	81-3138
	6-Hole Boxcar Plate	9 x 18 mm	1	81-3035
	8-Hole Boxcar Plate	11 x 25 mm	1	81-3036
	Y Plate	19 mm	1	81-3029
	X Plate	18 mm	1	81-3030
	T Plate	16 mm	1	81-3031
	Flexible Burr Hole Cover, Small	12 mm	1	81-2046
	Flexible Burr Hole Cover, Large	18 mm	1	81-2047
	Rigid Burr Hole Cover, Large	18 mm	1	81-2049
	Access Burr Hole Cover	18 mm	1	81-3181

Contact Bioplate to learn more about our complex and customizable plates not shown in this catalogue. 888-BIOPATE

1.9MM SCREWS




SELF-DRILLING SCREWS



Size	Qty	Code
3 mm long	1	81-2068
	5	81-2075
4 mm long	1	81-2069
	5	81-2076
5 mm long	1	81-2070
	5	81-2077
6 mm long	5	81-2078


1.9MM INSTRUMENTATION AND ACCESSORIES

SCREWDRIVER BLADE



Description	Qty	Code
Screwdriver Blade	1	81-3103

SCREWDRIVER HANDLE



Description	Qty	Code
Screwdriver Handle	1	81-3111

FORCEPS




Description	Qty	Code
Forceps - Angled	1	81-2024

CONTAINER



Description	Qty	Code
Plate & Screw Container	1	81-4206







1.9MM COMPACT PLUS STARTER KIT (81-3342)



Bioplate’s 1.9mm Compact Plus Starter Kit is a combination of the most frequently used plates, screws and instruments. It’s designed to be an introduction to our cranial fixation system. Delivered in our Compact Plus container, Bioplate provides a unique and flexible approach to rigid fixation. The Compact Plus container can be combined with other Bioplate containers for a stackable, compact footprint.

1.9MM STERILE KIT



Quantity Per Kit		Description	Code
3		2-hole Straight Plate, Large	81-3581 (BOX OF 10 KITS)
6		4 mm long Self-Drilling Screw	
Quantity Per Kit		Description	Code
1		X Plate	81-3586 (BOX OF 5 KITS)
6		4 mm long Self-Drilling Screw	
Quantity Per Kit		Description	Code
1		Flexible Burr Hole Cover, Large	81-3588 (BOX OF 5 KITS)
6		4 mm long Self-Drilling Screw	

Contact Bioplate to learn more about the contents of our **COMPACT PLUS KIT** and our customizable **STERILE KITS** not shown in this catalogue. 888-BIOPLATE

ORDER INFORMATION

888-BIOPATE

info@bioplate.com

PLATES - Ø1.9 mm System

81-3048	Straight Plate - 1.9 mm - 2 hole - 12 mm	Qty 1
81-3027	Straight Plate - 1.9 mm - 2 hole - 12 mm	Qty 5
81-3020	Bio-100 Plate - 1.9 mm - 2 hole - 15 mm	Qty 1
81-2037	Bio-100 Plate - 1.9 mm - 2 hole - 15 mm	Qty 5
81-3148	Straight Compression Plate - 1.9 mm - 2 hole - 15 mm	Qty 1
81-3028	Straight Plate - 1.9 mm - 4 hole - 20 mm	Qty 1
81-3034	Boxcar Plate - 1.9 mm - 4 hole - 11 x 11 mm	Qty 1
81-3138	Boxcar Plate - 1.9 mm - 4 hole - 15 x 15 mm	Qty 1
81-3035	Boxcar Plate - 1.9 mm - 6 hole - 9 x 18 mm	Qty 1
81-3036	Boxcar Plate - 1.9 mm - 8 hole - 11 x 25 mm	Qty 1
81-3029	Y Plate - 1.9 mm - 5 hole - 19 mm	Qty 1
81-3030	X Plate - 1.9 mm - 6 hole - 18 mm	Qty 1
81-3031	T Plate - 1.9 mm - 5 hole - 16 mm	Qty 1
81-3181	Access Burr Hole Cover - 1.9 mm - 5 hole - 18 mm	Qty 1
81-2046	Flexible Contoured Burr Hole Cover - 1.9 mm - 12 mm	Qty 1
81-2047	Flexible Contoured Burr Hole Cover - 1.9 mm - 18 mm	Qty 1
81-2049	Rigid Contoured Burr Hole Cover - 1.9 mm - 18 mm	Qty 1

STERILE KIT - Ø1.9 mm System

81-3581	Straight Plate - 1.9 mm - 2 hole - 15 mm 1.9 mm diam. Self-Drilling Screw - 4 mm long	(BOX OF 10 KITS)
81-3586	X Plate - 1.9 mm - 6 hole - 18 mm 1.9 mm diam. Self-Drilling Screw - 4 mm long	(BOX OF 5 KITS)
81-3588	Flexible Contoured Burr Hole Cover - 1.9 mm - 18 mm 1.9 mm diam. Self-Drilling Screw - 4 mm long	(BOX OF 5 KITS)

SCREWS - Ø1.9 mm System

81-2068	1.9 mm diam. Self-Drilling Screw - 3 mm long	Qty 1
81-2069	1.9 mm diam. Self-Drilling Screw - 4 mm long	Qty 1
81-2070	1.9 mm diam. Self-Drilling Screw - 5 mm long	Qty 1
81-2075	1.9 mm diam. Self-Drilling Screw - 3 mm long	Qty 5
81-2076	1.9 mm diam. Self-Drilling Screw - 4 mm long	Qty 5
81-2077	1.9 mm diam. Self-Drilling Screw - 5 mm long	Qty 5
81-2078	1.9 mm diam. Self-Drilling Screw - 6 mm long	Qty 5

ACCESSORIES - Ø1.9 mm System

81-3103	Screwdriver Blade - 1.9 mm diam.
81-3111	Screwdriver - Handle only
81-2024	Forceps - Angled
81-4206	1.9 mm Plate & Screw Container

