



BIOMESH®

CRANIAL FIXATION

This unique easy-to-use titanium mesh plate allows for precise, anatomic contouring to the cranial skeleton. Unlike standard sheeting, BioMesh® can be contoured into 3-dimensional shapes. This unique mesh can be molded, stretched or bent to fit any bony area. BioMesh® is available in an assortment of thicknesses and sizes for a variety of reconstructive needs.



BIOPLATE®


ADVANCING CRANIAL SURGERY

BIOMESH® PRODUCTS



BIOMESH®

	DESCRIPTION	QTY	CODE
2MC (0.45MM THICK)	BioMesh® 2MC - 33 x 48 mm	1	81-3021
	BioMesh® 2MC - 60 x 60 mm	1	81-3022
	BioMesh® 2MC - 82 x 97 mm	1	81-5788
	BioMesh® 2MC - 85 x 90 mm	1	81-3023
	BioMesh® 2MC - 132 x 134 mm	1	81-3140
3MC (0.6MM THICK)	BioMesh® 3MC - 36 x 52 mm	1	81-3024
	BioMesh® 3MC - 65 x 65 mm	1	81-3025
	BioMesh® 3MC - 90 x 92 mm	1	81-3026
	BioMesh® 3MC - 127 x 131 mm	1	81-3141
	BioMesh® 3MC - 200 x 200 mm	1	81-3142







SPECIALTY BIOMESH®

	DESCRIPTION	QTY	CODE
PTERIONAL PLATE	BioMesh® 2ME - 69 x 79 mm	1	81-3159
	BioMesh® Gap Plate - 29 x 88 mm	1	81-3184
BIOMESH® GAP PLATE	BioMesh® Gap Plate - 43 x 109 mm	1	81-3379
	4MC RECONSTRUCTION PLATE (0.8MM THICK)	BioMesh® 4MC - 200 x 200 mm	1



FEATURING HIGH-PRECISION TECHNOLOGY DESIGNED SPECIFICALLY FOR NEUROSURGERY, BIOPLATE® OFFERS A COMPREHENSIVE PRODUCT LINE AND SUPERIOR SERVICE TO THE SURGICAL COMMUNITY.

BIOMESH® ACCESSORIES

	DESCRIPTION	QTY	CODE
	BioMesh® Bender	1	81-3568
	Mesh Cutting Scissors	1	81-3380
	BioMesh® Container	1	81-4208 (1.5 mm)
		1	81-4209 (1.9 mm)
	Forceps - Angled	1	81-2024



BIOPLATE®

ADVANCING CRANIAL SURGERY

570 S Melrose St.
Placentia, CA 92870
info@bioplate.com
888-246-7528 (888-BIOPLATE)
bioplate.com



ORDER INFORMATION
888-BIOPLATE
WWW.BIOPLATE.COM

Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, warnings, and precautions.
Caution: US law restricts this device to sale by or on the order of a physician.
Bioplate® is a registered trademark of Bioplate, Inc. D00312 Rev D