




OSTEOPORE® BIORESORBABLE BONE SCAFFOLD

RESTORE BONE POST CRANIOTOMY




BIOPLATE®
ADVANCING CRANIAL SURGERY

BURR HOLE COVERS




OSTEOPORE® BURR HOLE COVER PC 22

CATALOG NO.	A, LENGTH (MM)	B, BREADTH (MM)	C, THICKNESS (MM)
PC22(11, 11, 4)	11	11	4
PC22(13, 13, 4)	13	13	4
PC22(14, 14, 4)	14	14	4
PC22(15, 15, 4)	15	15	4



OSTEOPORE® BURR HOLE COVER PC 23

CATALOG NO.	A, LENGTH (MM)	B, BREADTH (MM)	C, THICKNESS (MM)
PC23(11, 11, 4)	11	11	4
PC23(13, 13, 4)	13	13	4
PC23(14, 14, 4)	14	14	4
PC23(15, 15, 4)	15	15	4





OSTEOPORE® BURR HOLE COVER PC 24

CATALOG NO.	A, LENGTH (MM)	B, BREADTH (MM)	C, THICKNESS (MM)
PC24(11, 11, 4)	11	11	4
PC24(13, 13, 4)	13	13	4
PC24(14, 14, 4)	14	14	4
PC24(15, 15, 4)	15	15	4

STRIPS AND MESH

FILLING IN CRANIOTOMY GAPS CAN RESULT IN IMPROVED COSMESIS, REDUCED PAIN, AND INCREASED PATIENT SATISFACTION. THE OSTEOPORE IMPLANT FULLY RESORBS OVER 24 MONTHS. ^{1, 2, 3}



OSTEOPORE® STRIP PC 17

CATALOG NO.	A, LENGTH (MM)	B, BREADTH (MM)	C, THICKNESS (MM)
PC17(100, 2, 4)	100	2	4
PC17(100, 3.5, 4)	100	3.5	4

OSTEOPORE® MESH PC 11

CATALOG NO.	A, LENGTH (MM)	B, BREADTH (MM)	C, THICKNESS (MM)
PC11(50, 50, 1)	50	50	1
PC11(50, 50, 1.25)	50	50	1.25
PC11(50, 50, 2)	50	50	2



BIOPLATE®

ADVANCING CRANIAL SURGERY

570 S Melrose St.
Placentia, CA 92870
info@bioplate.com
888-246-7528 (888-BIOPLATE)
bioplate.com

Distributed by:
Bioplate, Inc.
Placentia, CA 92870

Manufactured by:
Osteopore International
Singapore 618305



ORDER INFORMATION
888-BIOPLATE
WWW.BIOPLATE.COM

¹Low SW et al. Use of Osteoplug polycaprolactone implants as novel burr-hole covers. Singapore Med J 2009; 50(8): 777-780.

²Schantz JT et al. Cranioplasty after trephination using a novel biodegradable burr hole cover: technical case report. Neurosurgery 58[ONS Suppl 1]:ONS-176, 2006.

³Yang M et al. Cranial reconstruction using a polycaprolactone implant after burr hole trephination. Journal of 3D Printing Medicine 2020; 4:9-16.

Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, warnings, and precautions.

Caution: US law restricts this device to sale by or on the order of a physician.

Osteopore, Osteopore International, and the stylized Osteopore logo are trademarks of Osteopore International Pte Ltd.

Bioplate is a registered trademark of Bioplate, Inc. D00551RevB